


The Combined Effect of Tibiofemoral Conformity & Tibial Base Design on Initial Stability of Cementless TKA



Authors: Chadd Clary¹, Hayden Wilson¹, James Deacy¹, Paul Rullkoetter¹
 Affiliations: ¹Center for Orthopaedic Biomechanics, University of Denver, Denver, CO


1. Introduction

- Initial fixation of cementless tibial baseplates is critical to long-term bony ingrowth
 - Key factors: bone quality, surgical technique, baseplate fixation features, and loading of the tray via the articulating surfaces
- Increased sagittal conformity between the insert and femur improves knee stability¹, but effects on tibial micromotion is unclear (Fig. 1)
- Goal: characterize the affect of tibiofemoral conformity, knee kinematics, and fixation features on tibial micromotion (Fig. 2)




3. Results

- Largest tibial micromotions were observed during stair descent and gait, following a consistent pattern.
- Increased micromotion correlated with lower conformity and higher condylar translations.




4. Discussion

- Experimental evidence suggests that articulation conformity is major driver of micromotion, but what about tray design factors?
- In FE model, swapped articulating surfaces:
 - For ATTUNE[®] Knee, swapping with a low conformity design increased micromotion over 30% for gait and stair descent⁵
 - For Competitor T, swapping with a high conformity design reduced micromotion between 30% - 40%⁵




2. Methods

- Tibial Sawbones[™] were implanted with three contemporary cementless tibial tray designs
 - Loaded via their respective femurs and inserts in an AMTI VIVO[™]
 - Gait, stair descent, and deep knee bend cycles²⁻⁴
 - Micromotions between the tray and bone were measured using Digital Image Correlation
- A finite-element (FE) model of the experiment was validated using results from the testing²
 - Isolated effects of tray fixation features from articular conformity by virtually exchanging the articulating features of the highest and lowest conformity designs and predicting the resulting micromotion.



3. Results (Cont.)

- Lowest conformity design had highest condylar translations during Gait and Stair Descent
 - Gait: 6.0mm (Attune) to 11.7mm (Competitor T)
 - Stair Descent : 6.6mm (Attune) compared to 11.7mm (Competitor T)
- Lowest conformity design had highest micromotions during Gait and Stair Descent
 - Gait: 101.4µm (Attune) to 146.5µm (Competitor T)
 - Stair Descent : 126.7µm (Attune) compared to 160.5µm (Competitor T)



5. Conclusion

- Limitations:
 - Experimental model did not have ligamentous or muscular constraint, which may improve stability (but loading does)
 - Foam bone has different mechanical properties than real bone
 - Initial fixation is a good predictor of long-term survivorship, but clinical studies and RSA evaluation still necessary to understand clinical performance
- Conclusion: Cementless tibial tray designs with optimized fixation features and moderately increased levels of tibial-femoral conformity can reduce A-P translations during activities of daily living, thus reducing micromotion and improving initial stability of cementless tibial trays.

References:
 1. Deacy J, Clary C, Wilson H, Rullkoetter P. The influence of tibial tray articular geometry on tibial motion stability: an experimental and finite element study. J Biomech. 2018 Apr 30;74(12):3151-7.
 2. Clary C, Wilson H, Deacy J, Rullkoetter P. Validation of model predicted tibial tray synthetic bone relative motion to cementless total knee arthroplasty during activities of daily living. J Biomech. 2018 Aug 15;77:115-124.
 3. Wilson H, Behrman S, Noyes A, Wang A, Rullkoetter P, Clary C. Tibial Micromotion during Activities of Daily Living. Orthopaedic Research Society 2018 Annual Meeting, New Orleans, LA, Poster 0775.
 4. Wilson H, Clary C, Deacy J, Rullkoetter P, Clary C. Tibial Base Micromotion for Active Posture and Locomotion during Activities of Daily Living. Biomechanics. 2020;55(1):1-10.
 5. Wilson H, Clary C, Deacy J, Rullkoetter P, Clary C. Tibial Base Micromotion for Active Posture and Locomotion during Activities of Daily Living. 2019 Annual Meeting, Austin, TX, Poster 1750.

183570-210719 DSUS
 Please refer to the instructions for use for a complete list of indications, contraindications, warnings, and precautions.
 ©DePuy Synthes 2021. All rights reserved.

The third-party trademarks used herein are the trademarks of their respective owners
 All products may not be available and/or approved or cleared by all global regulatory authorities.
 Please contact your sales representative for questions regarding regional product availability.

