



Reduce CRBSI: Always leave enough space for BIOPATCH® Protective Disk with CHG to be effective.



Secure catheter at least 0.5" (1.25 cm) from insertion site. This allows for proper placement of BIOPATCH Disk, which must have complete contact with skin for maximum efficacy.

Correct positioning of catheter also allows sufficient room for optimal cleaning of insertion site and under catheter hub to help minimize infections in patients.

BIOPATCH Disk is **PROVEN to reduce CRBSI in patients with central venous and arterial catheters.¹ Always make sure it is applied correctly.**

Don't secure catheter too close to entry point. This will prevent proper placement of BIOPATCH Disk and leave insufficient room to cleanse.



BIOPATCH Disk is **PROVEN** to reduce CRBSI.¹

- BIOPATCH Disk is the **ONLY** IV dressing with CHG proven in multiple, randomized controlled trials to reduce the incidence of catheter-related bloodstream infections (CRBSIs).
- **PROVEN** to reduce the incidence of CRBSI by 60%, and local infections by 44% in patients with central venous or arterial lines.¹
- Correct positioning of the catheter allows for cleaning under the catheter hub.

How to apply BIOPATCH Disk correctly:

1. Prepare insertion site according to hospital protocol.
2. Place BIOPATCH Disk around catheter **printed side up**.
3. Align radial slit with catheter (for easy removal).
4. Ensure slit edges touch (for maximum efficacy).
5. Secure catheter at least 0.5" (1.25 cm) from insertion site.
6. Apply transparent film dressing according to hospital protocol.



**For Full Prescribing Information or technical support
call 1-877-ETHICON (877-384-4266) or visit www.BIOPATCH.com**

ETHICON[®]

BIOPATCH[®] Products

Protect All Lines. Protect All Lives.[™]

*WARNING: Not for use on premature infants or patients with known sensitivity to CHG. Safety and effectiveness in children under 16 years of age has not been established.

References

1 Maki DG, Mermel L, Genthner D, Hua S, Chiacchierini RP. An evaluation of BIOPATCH[®] Antimicrobial Dressing compared to routine standard of care in the prevention of catheter-related bloodstream infection. Johnson & Johnson Wound Management, a division of ETHICON, INC. 2000. Data on file.