

Small Battery Drive II

Powerful, Durable and Comprehensive Battery Power Tool System



Small Battery Drive II

Battery power tool system for use in traumatology, hand and foot surgery

Powerful handpiece (85 Watts)
with highly durable motor

Simple mode selection from
forward to reverse and oscillation

Coupling for
quick and secure
attachment
exchange

Simple trigger
mechanisms for
easy speed
regulation and
mode control

14.4 V environment friendly*
Li-Ion battery offering high
battery capacity



Standing battery casing for
greater convenience

*As per Battery Directive 2006 / 66 / EC.

Rotating
attachments

K-Wire attachment



Extensive Portfolio of Attachments for a Wide Range of Applications

Added Power with the 14.4 V Li-Ion Battery Pack

85 Watts
14.4 V Li-Ion Battery Pack
532.103 / 532.132 / 532.104



Saw attachments

Ream attachments

Drill and screw
attachments

Saw attachment
for large crescentic
osteotomies

Innovative Reprocessing Solutions with designated washing basket





Universal Battery Charger II

The Universal Battery Charger II can charge and maintain all DePuy Synthes batteries and the TRS Power Module.



Washing and Sterilization Solution

Dedicated Washing Basket for optimal cleaning and sterilization conditions.



Service Solution

Our dedication to quality begins at the production of our valuable products but does not end there. We offer service contracts and customized solutions for each user. Annual maintenance of your Small Battery Drive II system is important to maintain safety, reliability and service life of your product.

Please also refer to the package insert(s) or other labeling associated with the devices identified in this brochure for additional information.

CAUTION: Federal Law restricts these devices to sale by or on the order of a physician.

Some devices listed in this brochure may not have been licensed in accordance with Canadian law and may not be for sale in Canada.

Please contact your sales consultant for items approved for sale in Canada.

Not all products may currently be available in all markets.



PART OF THE **Johnson & Johnson** FAMILY OF COMPANIES

Manufactured by:
Synthes USA Products, LLC
1302 Wrights Lane East
West Chester, PA 19380
USA

To order (USA): 800-327-6887
To order (Canada): 800-872-3080

Note: For recognized manufacturer, refer to the product label.

www.jnjmedicaldevices.com

© DePuy Synthes 2020. All rights reserved.
101032-200224 DSUS