

# Controlled compression for a strong anastomosis

## ETHICON™ Circular Stapler



### Designed with better healing in mind

- The ETHICON™ Circular Stapler with its Controlled Tissue Compression and delivered less tissue damage compared to fixed compression staplers.<sup>3</sup>

### More device feedback when firing

- Breakaway washer created more audible and tactile feedback than the Medtronic DST Series EEA.<sup>4</sup>

### Applicable specialties

- Colorectal
- Gastric
- Thoracic



Colon resection



Gastrectomy  
Gastric Bypass



Esophagectomy

### Less tissue displacement during closure and tissue strain during opening

- The ETHICON™ Circular Stapler with its continuously smooth anvil shaft created 88% less tissue movement during anvil closure compared to the Medtronic DST Series™ EEA™.<sup>1</sup>
- Fixed anvil design produced 73% less tissue strain during opening than the Medtronic DST Series™ EEA™ with its tilting anvil.<sup>2</sup>

## ETHICON™ Circular Stapler, XL Sealed



CODE	DIAMETER	# OF STAPLES	OPEN STAPLE LEG LENGTH	CLOSED STAPLE HEIGHT	SHAFT LENGTH	QUANTITY PER BOX
ECS21B	21mm	16	5.2mm	1.5 - 2.2mm	XL	3
ECS25B	25mm	20	5.2mm	1.5 - 2.2mm	XL	3
ECS29B	29mm	24	5.2mm	1.5 - 2.2mm	XL	3
ECS33B	33mm	28	5.2mm	1.5 - 2.2mm	XL	3

### How to order

All purchase orders are made to Johnson & Johnson Health Care Systems, Inc. (JJHCS) or your authorized distributor.

#### Electronic ordering options

(Note: Placing orders electronically avoids minimum order fees for your hospital.)

- **Johnson & Johnson Gateway: [www.jnjgateway.com/commerce](http://www.jnjgateway.com/commerce)** For questions about your order, visit the website or call 1-866-JNJ-GATE (565-4283)
- **Global Healthcare Exchange (GHX): [www.ghx.com](http://www.ghx.com)**  
For questions about your order, visit the website or call 1-800-YOUR-GHX (986-7449)
- **Electronic Data Interchange (EDI): JJHCS EDI**
- **Help Line: 1-800-262-2888**

#### Nonelectronic/manual ordering options

Call JJHCS at 1-800-255-2500 between 8:30 a.m. and 9:00 p.m. Eastern Standard Time or fax your order to 1-732-562-2212.

### Customer support

For product use assistance, clinical guidelines, service and repair, emergency assistance, copy of a 510(k) clearance letter or complaints, please contact our Customer Support Center at [customersupport@eesus.jnj.com](mailto:customersupport@eesus.jnj.com) or by calling **1-877-ETHICON** (384-4266).

Our Support Center is staffed 24 hours a day, 7 days a week by qualified nurses to answer your product-related questions.