

## STRATAFIX™ Knotless Tissue Control Devices

### Needles and Trocar Chart

Needle	Curvature	Length (mm)	Trocar Diameter (mm)				
			5	8	11	12	15
<b>RB-1</b>	1/2	15, 17, 30 16 x 16	No	Yes	Yes	Yes	Yes
<b>CT-3</b>	1/2	22 14 x 14	No	Yes	Yes	Yes	Yes
<b>SH-1</b>	1/2	22, 30	No	Yes	Yes	Yes	Yes
<b>CT-2</b>	1/2	15, 26, 30, 45	No	No	Yes	Yes	Yes
<b>CP-2</b>	1/2	26, 30, 45 14 x 14, 24 x 24	No	No	Yes	Yes	Yes
<b>SH</b>	1/2	15, 20, 26, 30, 45 10 x 10, 14 x 14	No	No	Yes	Yes	Yes
<b>OS-6</b>	1/2	36, 45 24 x 24	No	No	No	Yes	Yes
<b>MH</b>	1/2	36 7 x 7, 14 x 14, 24 x 24	No	No	No	Yes	Yes
<b>CT-1</b>	1/2	15, 26, 30, 45 14 x 14, 24 x 24	No	No	No	Yes	Yes
<b>CTB-1</b>	1/2	36, 45	No	No	No	Yes	Yes
<b>MO-4</b>	1/2	36	No	No	No	No	Yes
<b>CT</b>	1/2	40	No	No	No	No	Yes
<b>OS-8</b>	1/2	40	No	No	No	No	Yes
<b>CTX</b>	1/2	48	No	No	No	No	No
<b>CTXB</b>	1/2	48	No	No	No	No	No
<b>V-7</b>	1/2	30	No	Yes	Yes	Yes	Yes
<b>V-34</b>	1/2	30	No	No	No	Yes	Yes
<b>V-4</b>	3/8	17 24 x 24	No	Yes	Yes	Yes	Yes
<b>PS-2</b>	3/8	19	No	Yes	Yes	Yes	Yes
<b>FS-2</b>	3/8	19	No	Yes	Yes	Yes	Yes
<b>PS-1</b>	3/8	24	No	Yes	Yes	Yes	Yes
<b>FS-1</b>	3/8	24	No	Yes	Yes	Yes	Yes
<b>FS</b>	3/8	26	No	Yes	Yes	Yes	Yes
<b>V-26</b>	3/8	14 x 14, 24 x 24	No	Yes	Yes	Yes	Yes
<b>PSL</b>	3/8	30	No	No	Yes	Yes	Yes
<b>BB</b>	3/8	15, 30	No	Yes	Yes	Yes	Yes
<b>ST-4</b>	Straight	15, 30	Yes	Yes	Yes	Yes	Yes
<b>EN-S</b>	Half-curved	30	Yes	Yes	Yes	Yes	Yes
<b>UR-6</b>	5/8	30	No	No	Yes	Yes	Yes