

Inventory Control Form

3.0 mm, 4.0 mm Titanium Cannulated Screws



Patient Information:

Synthes (USA)
To Order: (800) 523-0322

Synthes (Canada) Ltd.
To Order: (800) 668-1119

Date: _____

Hospital: _____

Surgeon: _____

Procedure: _____

Implants

3.0 mm Titanium Cannulated Screws, short thread



	LENGTH		LENGTH
<input type="checkbox"/>	402.608 8 mm	<input type="checkbox"/>	402.622 22 mm
<input type="checkbox"/>	402.609 9 mm	<input type="checkbox"/>	402.623 23 mm
<input type="checkbox"/>	402.610 10 mm	<input type="checkbox"/>	402.624 24 mm
<input type="checkbox"/>	402.611 11 mm	<input type="checkbox"/>	402.625 25 mm
<input type="checkbox"/>	402.612 12 mm	<input type="checkbox"/>	402.626 26 mm
<input type="checkbox"/>	402.613 13 mm	<input type="checkbox"/>	402.627 27 mm
<input type="checkbox"/>	402.614 14 mm	<input type="checkbox"/>	402.628 28 mm
<input type="checkbox"/>	402.615 15 mm	<input type="checkbox"/>	402.629 29 mm
<input type="checkbox"/>	402.616 16 mm	<input type="checkbox"/>	402.630 30 mm
<input type="checkbox"/>	402.617 17 mm	<input type="checkbox"/>	402.632 32 mm
<input type="checkbox"/>	402.618 18 mm	<input type="checkbox"/>	402.634 34 mm
<input type="checkbox"/>	402.619 19 mm	<input type="checkbox"/>	402.636 36 mm
<input type="checkbox"/>	402.620 20 mm	<input type="checkbox"/>	402.638 38 mm
<input type="checkbox"/>	402.621 21 mm	<input type="checkbox"/>	402.640 40 mm

4.0 mm Titanium Cannulated Screws, short thread



	LENGTH		LENGTH
<input type="checkbox"/>	407.610 10 mm	<input type="checkbox"/>	407.640 40 mm
<input type="checkbox"/>	407.612 12 mm	<input type="checkbox"/>	407.642 42 mm
<input type="checkbox"/>	407.614 14 mm	<input type="checkbox"/>	407.644 44 mm
<input type="checkbox"/>	407.616 16 mm	<input type="checkbox"/>	407.646 46 mm
<input type="checkbox"/>	407.618 18 mm	<input type="checkbox"/>	407.648 48 mm
<input type="checkbox"/>	407.620 20 mm	<input type="checkbox"/>	407.650 50 mm
<input type="checkbox"/>	407.622 22 mm	<input type="checkbox"/>	407.652 52 mm
<input type="checkbox"/>	407.624 24 mm	<input type="checkbox"/>	407.654 54 mm
<input type="checkbox"/>	407.626 26 mm	<input type="checkbox"/>	407.656 56 mm
<input type="checkbox"/>	407.628 28 mm	<input type="checkbox"/>	407.658 58 mm
<input type="checkbox"/>	407.630 30 mm	<input type="checkbox"/>	407.660 60 mm
<input type="checkbox"/>	407.632 32 mm	<input type="checkbox"/>	407.664 64 mm
<input type="checkbox"/>	407.634 34 mm	<input type="checkbox"/>	407.668 68 mm
<input type="checkbox"/>	407.636 36 mm	<input type="checkbox"/>	407.672 72 mm
<input type="checkbox"/>	407.638 38 mm		

3.0 mm Titanium Cannulated Screws, long thread



	LENGTH		LENGTH
<input type="checkbox"/>	402.714 14 mm	<input type="checkbox"/>	402.725 25 mm
<input type="checkbox"/>	402.715 15 mm	<input type="checkbox"/>	402.726 26 mm
<input type="checkbox"/>	402.716 16 mm	<input type="checkbox"/>	402.727 27 mm
<input type="checkbox"/>	402.717 17 mm	<input type="checkbox"/>	402.728 28 mm
<input type="checkbox"/>	402.718 18 mm	<input type="checkbox"/>	402.729 29 mm
<input type="checkbox"/>	402.719 19 mm	<input type="checkbox"/>	402.730 30 mm
<input type="checkbox"/>	402.720 20 mm	<input type="checkbox"/>	402.732 32 mm
<input type="checkbox"/>	402.721 21 mm	<input type="checkbox"/>	402.734 34 mm
<input type="checkbox"/>	402.722 22 mm	<input type="checkbox"/>	402.736 36 mm
<input type="checkbox"/>	402.723 23 mm	<input type="checkbox"/>	402.738 38 mm
<input type="checkbox"/>	402.724 24 mm	<input type="checkbox"/>	402.740 40 mm

4.0 mm Titanium Cannulated Screws, long thread



	LENGTH		LENGTH
<input type="checkbox"/>	407.716 16 mm	<input type="checkbox"/>	407.742 42 mm
<input type="checkbox"/>	407.718 18 mm	<input type="checkbox"/>	407.744 44 mm
<input type="checkbox"/>	407.720 20 mm	<input type="checkbox"/>	407.746 46 mm
<input type="checkbox"/>	407.722 22 mm	<input type="checkbox"/>	407.748 48 mm
<input type="checkbox"/>	407.724 24 mm	<input type="checkbox"/>	407.750 50 mm
<input type="checkbox"/>	407.726 26 mm	<input type="checkbox"/>	407.752 52 mm
<input type="checkbox"/>	407.728 28 mm	<input type="checkbox"/>	407.754 54 mm
<input type="checkbox"/>	407.730 30 mm	<input type="checkbox"/>	407.756 56 mm
<input type="checkbox"/>	407.732 32 mm	<input type="checkbox"/>	407.758 58 mm
<input type="checkbox"/>	407.734 34 mm	<input type="checkbox"/>	407.760 60 mm
<input type="checkbox"/>	407.736 36 mm	<input type="checkbox"/>	407.764 64 mm
<input type="checkbox"/>	407.738 38 mm	<input type="checkbox"/>	407.768 68 mm
<input type="checkbox"/>	407.740 40 mm	<input type="checkbox"/>	407.772 72 mm

Inventory Control Form

3.0 mm, 4.0 mm Titanium Cannulated Screws

Implants (cont'd)

Washers for 3.0 mm Titanium Cannulated Screws

419.89 Titanium Threaded Washer



419.972 Titanium Washer, 6.5 mm



Washers for 4.0 mm Titanium Cannulated Screws

419.98 Titanium Washer, 7.0 mm



Instruments

Instruments for 3.0 mm Titanium Cannulated Screws

- 292.622 1.1 mm Threaded Guide Wire, 150 mm
- 292.623 1.1 mm Non-Threaded Guide Wire, 150 mm
- 310.221 2.0 mm Cannulated Drill Bit, quick coupling, 150 mm

Instruments for 4.0 mm Titanium Cannulated Screws

- 310.67 2.7 mm Cannulated Drill Bit, quick coupling, 160 mm
- 311.63 Cannulated Tap for 4.0 mm Cannulated Screws, 147 mm
- 900.722 1.25 mm Threaded Guide Wire, 150 mm
- 900.721 1.25 mm Non-Threaded Guide Wire, 150 mm