


TWISTR™ RETROGRADE REAMER & CRUCIATE+ INSTRUMENT SYSTEM

‘One-Size-Fits-All’ Reaming Combined with Precision Aiming¹



The TWISTR™ Retrograde Reamer blade can be adjusted to drill bone tunnels between 6mm & 12mm, including half sizes

Meaningful Differentiation

- TWISTR Retrograde Reamer: One innovative device capable of drilling 13 different tunnel sizes
 - One-size-fits-all design eliminates the need to stock multiple reamer sizes, streamlining inventory management
 - Can be adjusted mid-case to drill multiple tunnel sizes, increasing intra-operative versatility and reducing cost
- MITEK Cruciate+ Instrument System increases drilling accuracy for precision, anatomic tunnel placement¹





Getting your patients back to their passion is our passion.



Ordering Information

	Product Code	Description
Disposables	232000	TWISTR™ Retrograde Reamer Single Patient Use Disposable Retrograde Reamer, 6-12mm
	232301	Outside-In ACL Disposables Kit Disposable instrumentation kit for outside-in drilling during ACL reconstruction. Kit contains: outside-in suture passing pin, 1.1mm guidewire, 14" drill-tip style beath pin, malleable graft retractor, marking pen, and ruler
	Product Code	Description
Tray	215203	MITEK Cruciate+ Instrument System Tray Sterilization Tray
	215205	MITEK Cruciate+ Instrument System Tray Lid Sterilization Tray Lid

	Product Code	Description
Carriage, Bullets & Inserts	219001	Guide Carriage Side & Back Loading Guide Carriage
	219002	Obturator for Retrograde Reamer Bullet 4.8mm Obturator for Use with 219004 Bullet
	219003	Beath Pin Bullet 2.4mm Inner Diameter Bullet, for use with Standard Beath Pin
	219004	Retrograde Reamer Bullet 4.8mm Inner Diameter Bullet, for use with TWISTR Reamer
	219005	2.4mm Insert for Retrograde Reamer Bullet 2.4mm Beath Pin Insert / Adapter for 219004 Bullet
	219018	Outside-In Suture Passing Pin Forked Tip Reusable Suture Passer

	Tibial		Femoral	
Product Code	219006	219007	219008	219009
Description	ACL Tibial Point Aimer Point-to-Point Style Tibial Aimer	ACL Tibial Capture Aimer Enclosed (Capture-style) Tibial Aimer	ACL/PCL Femoral Point Aimer Point-to-Point Style Tibial Aimer	ACL/PCL Femoral Capture Aimer Enclosed (Capture-style) Femoral Aimer
Aimer Arms				
				

Limited Warranty and Disclaimer: DePuy Mitek, Inc. products are sold with a limited warranty to the original purchaser against defects in workmanship and materials. Any other express or implied warranties, including warranties of merchantability or fitness, are hereby disclaimed.

WARNING: In the USA, this product has labeling limitations. See package insert for complete information.

CAUTION: USA Law restricts these devices to sale by or on the order of a physician.

The third party trademarks used herein are the trademarks of their respective owners.



Medos International SÄRL
Chemin-Blanc 38
2400 Le Locle
Switzerland

www.depuyssynthes.com

Distributed in the USA by:
DePuy Mitek, Inc.
325 Paramount Drive
Raynham, MA 02767 USA
Tel: +1 (800) 382-4682

