

Inventory Control Form

Distal Radius Sterile Kit

Case Information (eg, Case ID):

DO NOT PUT ANY
PATIENT IDENTIFIABLE
INFORMATION IN
THIS FIELD

To order (USA): 800-523-0322

To order (Canada): 855-946-8999

Date: _____

Hospital: _____

Surgeon: _____

Procedure: _____

Core Kits

Titanium

2.4 mm VA LCP® Two-Column Volar Distal Radius Plates, Implant and Instrument Kit, sterile*

	Head size (mm)	Head holes	Shaft holes	Length (mm)	Left/Right
<input type="checkbox"/> 04.111.530KS	narrow 19.5	6	3	51	R
<input type="checkbox"/> 04.111.531KS	narrow 19.5	6	3	51	L
<input type="checkbox"/> 04.111.630KS	standard 22	6	3	54	R
<input type="checkbox"/> 04.111.631KS	standard 22	6	3	54	L
<input type="checkbox"/> 04.111.730KS	wide 25.5	7	3	55	R
<input type="checkbox"/> 04.111.731KS	wide 25.5	7	3	55	L

Stainless Steel

2.4 mm VA LCP Two-Column Volar Distal Radius Plates, Implant and Instrument Kit, sterile**

	Head size (mm)	Head holes	Shaft holes	Length (mm)	Left/Right
<input type="checkbox"/> 02.111.530KS	narrow 19.5	6	3	51	R
<input type="checkbox"/> 02.111.531KS	narrow 19.5	6	3	51	L
<input type="checkbox"/> 02.111.630KS	standard 22	6	3	54	R
<input type="checkbox"/> 02.111.631KS	standard 22	6	3	54	L
<input type="checkbox"/> 02.111.730KS	wide 25.5	7	3	55	R
<input type="checkbox"/> 02.111.731KS	wide 25.5	7 </td <td>3</td> <td>55</td> <td>L</td>	3	55	L



*Plates are made from commercially pure titanium (TiCP), and screws are made from titanium alloy (Ti-6Al-7Nb).

**Plates and screws are made from stainless steel (316L).

Some devices in this form may not have been licensed in accordance with Canadian law and may not be for sale in Canada. Please contact your sales consultant for items approved for sale in Canada.


Inventory Control Form

Distal Radius Sterile Kit

Instruments

Instrument-Only Kit


03.111.100S 2.4 mm VA LCP Instrument Kit, sterile



Templates

Sizing Templates

03.111.120S 2.4 mm VA LCP 2-column Volar Distal Radius Plate Templates, sterile



03.111.117S Implant Removal Screwdriver



Some devices in this form may not have been licensed in accordance with Canadian law and may not be for sale in Canada. Please contact your sales consultant for items approved for sale in Canada.



Manufactured or distributed by:

Synthes USA Products, LLC
1302 Wrights Lane East
West Chester, PA 19380

Synthes USA, LLC
1101 Synthes Avenue
Monument, CO 80132

To order (USA): 800-523-0322
To order (Canada): 855-946-8999

Note: For recognized manufacturer, refer to the product label.

www.depuysynthes.com