

ECHELON CIRCULAR™ Powered Stapler **Competitive Comparison**

Ethicon ECHELON CIRCULAR™ Powered Stapler

Medtronic DST Series™ EEA™



	CDH23P	CDH25P	CDH29P	CDH31P	EEA21	EEA2135	EEA25	EEA2535	EEA28	EEA2835	EEA31	EEA33
DEVICE												
Size	23mm	25mm	29mm	31mm	21mm		25mm		28mm		31mm	33mm
Staple rows	2	2	2	2	2		2		2		2	2
Number of staples	20	22	26	30	18		22		26		30	32
Open staple height	5.2mm	5.2mm	5.2mm	5.2mm	4.8mm	3.5mm	4.8mm	3.5mm	4.8mm	3.5mm	4.8mm	3.5mm
Nominal closed staple height	1.5 - 2.2mm	1.5 - 2.2mm	1.5 - 2.2mm	1.5 - 2.2mm	2.0mm	1.5mm	2.0mm	1.5mm	2.0mm	1.5mm	2.0mm	2.0mm
Cutting diameter	14.6mm	16.5mm	20.4mm	22.5mm	12.5mm		16.6mm		19.5mm		22.5mm	24.6mm
Closure system	Controlled tissue compression Firing range window Variable height staples Dual-stage closing mechanism				Zero-gap anvil Ready to fire indicator window Fixed height staples Dual-stage closing mechanism							
Cutting washer	Rigid washer with audible & tactile indication of firing completion				Thin film							
Cartridge	Gripping Surface Technology				Standard flat deck							
Staple type	Titanium alloy 3D-form Round wire				Titanium B-form Rectangular wire with Directional Stapling Technology							

The information provided is based on the products' respective IFUs and websites, and Ethicon internal testing.

The third-party trademarks used herein are trademarks of their respective owners.

For complete product information, including full steps for use, indications, contraindications, warnings and precautions, please see the Instructions for Use.

©2017 Ethicon, Inc. All rights reserved. 071477-170419