

INVENTORY CONTROL FORM

2.4 mm LCP® Distal Radius Plate System



COMPANIES OF Johnson & Johnson

People inspired™

Case Information:

DO NOT PUT ANY
PATIENT IDENTIFIABLE
INFORMATION IN
THIS FIELD

Synthes USA, LLC

To order: (800) 523-0322

For Canada

To order: (800) 668-1119

Date: _____

Hospital: _____

Surgeon: _____

Procedure: _____

14-digit code: _____

Implants

2.4 mm Locking Screws, self-tapping, T8 STARDRIVE™ Screwdriver Recess




	LENGTH		LENGTH
<input type="checkbox"/>	212.806 6 mm	<input type="checkbox"/>	212.820 20 mm
<input type="checkbox"/>	212.808 8 mm	<input type="checkbox"/>	212.822 22 mm
<input type="checkbox"/>	212.810 10 mm	<input type="checkbox"/>	212.824 24 mm
<input type="checkbox"/>	212.812 12 mm	<input type="checkbox"/>	212.826 26 mm
<input type="checkbox"/>	212.814 14 mm	<input type="checkbox"/>	212.828 28 mm
<input type="checkbox"/>	212.816 16 mm	<input type="checkbox"/>	212.830 30 mm
<input type="checkbox"/>	212.818 18 mm		

2.4 mm Cortex Screws, self-tapping, T8 STARDRIVE Screwdriver Recess



	LENGTH		LENGTH
<input type="checkbox"/>	201.756 6 mm	<input type="checkbox"/>	201.770 20 mm
<input type="checkbox"/>	201.758 8 mm	<input type="checkbox"/>	201.772 22 mm
<input type="checkbox"/>	201.760 10 mm	<input type="checkbox"/>	201.774 24 mm
<input type="checkbox"/>	201.762 12 mm	<input type="checkbox"/>	201.776 26 mm
<input type="checkbox"/>	201.764 14 mm	<input type="checkbox"/>	201.778 28 mm
<input type="checkbox"/>	201.766 16 mm	<input type="checkbox"/>	201.780 30 mm
<input type="checkbox"/>	201.768 18 mm		

2.7 mm Cortex Screws, self-tapping, T8 STARDRIVE Screwdriver Recess



	LENGTH		LENGTH
<input type="checkbox"/>	202.870 10 mm	<input type="checkbox"/>	202.882 22 mm
<input type="checkbox"/>	202.872 12 mm	<input type="checkbox"/>	202.884 24 mm
<input type="checkbox"/>	202.874 14 mm	<input type="checkbox"/>	202.886 26 mm
<input type="checkbox"/>	202.876 16 mm	<input type="checkbox"/>	202.888 28 mm
<input type="checkbox"/>	202.878 18 mm	<input type="checkbox"/>	202.890 30 mm
<input type="checkbox"/>	202.880 20 mm		

1.8 mm LCP Buttress Pins



	LENGTH		LENGTH
<input type="checkbox"/>	200.190 12 mm	<input type="checkbox"/>	200.195 22 mm
<input type="checkbox"/>	200.191 14 mm	<input type="checkbox"/>	200.196 24 mm
<input type="checkbox"/>	200.192 16 mm	<input type="checkbox"/>	200.197 26 mm
<input type="checkbox"/>	200.193 18 mm	<input type="checkbox"/>	200.198 28 mm
<input type="checkbox"/>	200.194 20 mm	<input type="checkbox"/>	200.199 30 mm

2.4 mm LCP® Volar Distal Radius Plates, 5 holes head



		SHAFT HOLES	LENGTH
<input type="checkbox"/>	242.491 Left	3	43 mm
<input type="checkbox"/>	242.492 Left	5	61 mm
<input type="checkbox"/>	242.493 Right	3	43 mm
<input type="checkbox"/>	242.494 Right	5	61 mm

2.4 mm LCP Volar Distal Radius Plates, extra-articular, 5 holes head



		SHAFT HOLES	LENGTH
<input type="checkbox"/>	242.458 Right	3	48 mm
<input type="checkbox"/>	242.459 Right	5	66 mm
<input type="checkbox"/>	242.461 Left	3	48 mm
<input type="checkbox"/>	242.462 Left	5	66 mm

2.4 mm LCP Volar Distal Radius Plates, extra-articular, 4 holes head



		SHAFT HOLES	LENGTH
<input type="checkbox"/>	242.464 Right	3	47 mm
<input type="checkbox"/>	242.465 Right	5	65 mm
<input type="checkbox"/>	242.467 Left	3	47 mm
<input type="checkbox"/>	242.468 Left	5	65 mm

2.4 mm LCP Distal Radius Plates, straight



		SHAFT HOLES	LENGTH
<input type="checkbox"/>	242.479	5	48 mm
<input type="checkbox"/>	242.490	6	57 mm

INVENTORY CONTROL FORM

2.4 mm LCP Distal Radius Plate System



COMPANIES OF Johnson & Johnson

People inspired™

Implants (cont'd)

2.4 mm LCP Distal Radius L-Plates, 2 holes head



			SHAFT HOLES	LENGTH
<input type="checkbox"/>	242.500	-90°	3	40 mm
<input type="checkbox"/>	242.501	-90°	4	49 mm
<input type="checkbox"/>	242.502	+90°	3	40 mm
<input type="checkbox"/>	242.503	+90°	4	49 mm

2.4 mm LCP Distal Radius L-Plates, 3 holes head



			SHAFT HOLES	LENGTH
<input type="checkbox"/>	242.504	+90°	3	40 mm
<input type="checkbox"/>	242.505	+90°	4	49 mm
<input type="checkbox"/>	242.506	-90°	3	40 mm
<input type="checkbox"/>	242.507	-90°	4	49 mm

<input type="checkbox"/>	242.508	+20°	3	43 mm
<input type="checkbox"/>	242.509	+20°	4	52 mm
<input type="checkbox"/>	242.511	-20°	3	43 mm
<input type="checkbox"/>	242.512	-20°	4	52 mm



2.4 mm LCP Distal Radius T-Plates, 3 holes head



			SHAFT HOLES	LENGTH
<input type="checkbox"/>	242.477		3	40 mm
<input type="checkbox"/>	242.478		4	49 mm

2.4 mm LCP Volar Column Distal Radius Plates, 6 holes head[◇]



			SHAFT HOLES	LENGTH
<input type="checkbox"/>	02.110.730	Right	3	45 mm
<input type="checkbox"/>	02.110.731	Left	3	45 mm
<input type="checkbox"/>	02.110.740	Right	4	54 mm
<input type="checkbox"/>	02.110.741	Left	4	54 mm
<input type="checkbox"/>	02.110.750	Right	5	63 mm
<input type="checkbox"/>	02.110.751	Left	5	63 mm

2.4 mm LCP Volar Column Distal Radius Plates, 7 holes head[◇]



			SHAFT HOLES	LENGTH
<input type="checkbox"/>	02.110.830	Right	3	45 mm
<input type="checkbox"/>	02.110.831	Left	3	45 mm
<input type="checkbox"/>	02.110.840	Right	4	54 mm
<input type="checkbox"/>	02.110.841	Left	4	54 mm
<input type="checkbox"/>	02.110.850	Right	5	63 mm
<input type="checkbox"/>	02.110.851	Left	5	63 mm

2.4 mm LCP Volar Column Distal Radius Plates, 8 holes head



			SHAFT HOLES	LENGTH
<input type="checkbox"/>	02.110.430	Right	3	50.5 mm
<input type="checkbox"/>	02.110.431	Left	3	50.5 mm
<input type="checkbox"/>	02.110.440	Right	4	59.5 mm
<input type="checkbox"/>	02.110.441	Left	4	59.5 mm
<input type="checkbox"/>	02.110.450	Right	5	68.5 mm
<input type="checkbox"/>	02.110.451	Left	5	68.5 mm

2.4 mm LCP Volar Column Distal Radius Plates, 9 holes head



			SHAFT HOLES	LENGTH
<input type="checkbox"/>	02.110.330	Right	3	50.5 mm
<input type="checkbox"/>	02.110.331	Left	3	50.5 mm
<input type="checkbox"/>	02.110.340	Right	4	59.5 mm
<input type="checkbox"/>	02.110.341	Left	4	59.5 mm
<input type="checkbox"/>	02.110.350	Right	5	68.5 mm
<input type="checkbox"/>	02.110.351	Left	5	68.5 mm

LCP Dia-Meta Volar Distal Radius Plates[◇]



			SHAFT HOLES
<input type="checkbox"/>	02.110.005	Right	5
<input type="checkbox"/>	02.110.007	Right	7
<input type="checkbox"/>	02.110.009	Right	9
<input type="checkbox"/>	02.110.011	Right	11
<input type="checkbox"/>	02.110.105	Left	5
<input type="checkbox"/>	02.110.107	Left	7
<input type="checkbox"/>	02.110.109	Left	9
<input type="checkbox"/>	02.110.111	Left	11
<input type="checkbox"/>	02.110.013S	Right, sterile only	13
<input type="checkbox"/>	02.110.015S	Right, sterile only	15
<input type="checkbox"/>	02.110.113S	Left, sterile only	13
<input type="checkbox"/>	02.110.115S	Left, sterile only	15

Kirschner Wires, with trocar point, 150 mm

<input type="checkbox"/>	292.12	1.25 mm, 10/pkg.
<input type="checkbox"/>	292.16	1.6 mm, 10/pkg.
<input type="checkbox"/>	02.110.300	1.8 mm, 10/pkg.

[◇] Available nonsterile or sterile-packed. Add "S" to catalog number to order sterile product.

INVENTORY CONTROL FORM

2.4 mm LCP Distal Radius Plate System



COMPANIES OF Johnson & Johnson

People inspired™

Instruments

Drill Bits



<input type="checkbox"/>	310.19	2.0 mm, quick coupling, 100 mm
<input type="checkbox"/>	310.26	2.7 mm, quick coupling, 100 mm
<input type="checkbox"/>	310.509	1.8 mm with depth mark, quick coupling, 110 mm
<input type="checkbox"/>	310.530	2.4 mm, quick coupling, 100 mm

1.25 mm Plate Reduction Wires, threaded tip, 150 mm



02.111.500.10 Small stop



02.111.501.10 Large stop

Threaded Drill Guide

323.035 1.8 mm, short, 5/pkg.

