

INTRAFIX® ADVANCE
Tibial Fastener System



Legacy of Strength in Tibial Fixation



INTRAFIX® ADVANCE Tibial Fastener System offers trusted performance with simplicity and speed.

- Redesigned implants and instrumentation for greater ease of use and reliability
- High fixation strength and low displacement¹
- Sheath protects graft and provides 360° of graft-to-bone compression
- Available in BIOCRYL® RAPIDE™ biocomposite material and PEEK
- Simplified implant sizing scheme

Biomechanical performance that surgeons rely on



Product #	UOM	IMPLANTS BIOCRYL RAPIDE (Screw & Sheath)
254800	EA	INTRAFIX ADVANCE BR Screw, 7 x 23 mm & Small 23 mm Sheath
254801	EA	INTRAFIX ADVANCE BR Screw, 8 x 23 mm & Small 23 mm Sheath
254802	EA	INTRAFIX ADVANCE BR Screw, 9 x 23 mm & Large 23 mm Sheath
254803	EA	INTRAFIX ADVANCE BR Screw, 10 x 23 mm & Large 23 mm Sheath
254806	EA	INTRAFIX ADVANCE BR Screw, 7 x 30 mm & Small 30 mm Sheath
254807	EA	INTRAFIX ADVANCE BR Screw, 8 x 30 mm & Small 30 mm Sheath
254808	EA	INTRAFIX ADVANCE BR Screw, 9 x 30 mm & Large 30 mm Sheath
254809	EA	INTRAFIX ADVANCE BR Screw, 10 x 30 mm & Large 30 mm Sheath
IMPLANTS PEEK (Screw only)		
254814	EA	INTRAFIX ADVANCE PEEK Screw, 7 x 23 mm
254815	EA	INTRAFIX ADVANCE PEEK Screw, 8 x 23 mm
254816	EA	INTRAFIX ADVANCE PEEK Screw, 9 x 23 mm
254817	EA	INTRAFIX ADVANCE PEEK Screw, 10 x 23 mm
254820	EA	INTRAFIX ADVANCE PEEK Screw, 7 x 30 mm
254821	EA	INTRAFIX ADVANCE PEEK Screw, 8 x 30 mm
254822	EA	INTRAFIX ADVANCE PEEK Screw, 9 x 30 mm
254823	EA	INTRAFIX ADVANCE PEEK Screw, 10 x 30 mm
254824	EA	INTRAFIX ADVANCE PEEK Screw, 11 x 30 mm
254841	EA	INTRAFIX ADVANCE PEEK Screw, 12 x 30 mm
IMPLANTS POLYPROPYLENE (Sheath only)		
254818	EA	INTRAFIX ADVANCE PP Sheath, Small 23 mm
254819	EA	INTRAFIX ADVANCE PP Sheath, Large 23 mm
254825	EA	INTRAFIX ADVANCE PP Sheath, Small 30 mm
254826	EA	INTRAFIX ADVANCE PP Sheath, Large 30 mm
254827	EA	INTRAFIX ADVANCE PP Sheath, Extra Large 30 mm

Product #	UOM	INSTRUMENTS
254828	EA	Trial/Dilator - Small
254829	EA	Trial/Dilator - Large
254830	EA	Trial/Dilator - Extra Large
254831	EA	Sheath Inserter - Small 23 mm
254832	EA	Sheath Inserter - Large 23 mm
254833	EA	Sheath Inserter - Small 30 mm
254834	EA	Sheath Inserter - Large 30 mm
254835	EA	Sheath Inserter - Extra Large 30 mm
254836	EA	Hex Head Screw Driver
254847	EA	Ratchet Handle
254838	EA	Tensioner
254842	EA	Slap Hammer
254839	EA	Sterilization Case (Outer)
254845	EA	Sterilization (Inner) Tray 1
254846	EA	Sterilization (Inner) Tray 2
254848	EA	Sterilization Lid
254844	EA	Guidewire 12" Sterile
254843	PK	Guidewire 12" (6 pack) Sterile



For more information, contact your Mitek Sports Medicine sales representative, or call 1-800-382-4682.



PART OF THE **Johnson & Johnson** FAMILY OF COMPANIES

DePuy Synthes Mitek Sports Medicine
 325 Paramount Drive
 Raynham, MA 02767-0350
 T. +1 (800) 382-4682

www.depuyssynthes.com

The third-party trademarks used herein are the trademarks of their respective owners.
 Reference Internal Data: 1. DHF-103030262.

Limited Warranty and Disclaimer: DePuy Mitek, Inc. products are sold with a limited warranty to the original purchaser against defects in workmanship and materials. Any other express or implied warranties, including warranties of merchantability or fitness, are hereby disclaimed.

WARNING: In the USA, this product has labeling limitations. See package insert for complete information.

CAUTION: USA Law restricts these devices to sale by or on the order of a physician.

To order in the U.S. call Johnson & Johnson Health Care Systems Inc. Customer Support Services at 800-255-2500.

Not all products are currently available in all markets.

Medos International SÄRL
 Chemin-Blanc 38
 2400 Le Locle
 Switzerland

