

## Versatile and easy to use in a variety of electrosurgical procedures

MEGADYNE™ MEGA SOFT™ Reusable Patient Return Electrodes meet many OR needs

Product	MEGA SOFT™	MEGA SOFT™ Universal Plus	MEGA SOFT™ Universal
			
<b>Pad size</b>	46 inches long x 20 inches wide	36 inches long x 20 inches wide	36 inches long x 20 inches wide
<b>Pad thickness</b>	1/2-inch	1/4-inch	1/8-inch
<b>Key features</b>	Pressure reduction	Patient comfort, lightweight, pressure reduction	Versatile and lightweight
<b>Patient size</b>	25 lbs. and up	0.8 lbs. and up	0.8 lbs. and up
<b>Pad weight</b>	18.5 lbs.	8 lbs.	4 lbs.
<b>Pad use</b>	One-sided	One-sided	Two-sided
<b>Lightweight</b>		✓	✓
<b>Flexible</b>		✓	✓
<b>Pressure reduction</b>	✓	✓	
<b>Comfort</b>	✓	✓	
<b>Max life</b>	24 months	24 months	24 months
<b>Additional information</b>	Single cord option. Dual cord option for use with two generators.	Single cord option. Dual cord option for use with two generators. Pediatrics to adults.	Single cord option. Dual cord option for use with two generators. Pediatrics to adults.

Compatible with most major brands of isolated electrosurgical generators.