

Versatile and easy to use in a variety of electrosurgical procedures¹

MEGADYNE™ MEGA SOFT™ Reusable Patient Return Electrodes meet many OR needs

Product	MEGA SOFT™	MEGA SOFT™ Universal Plus	MEGA SOFT™ Universal
			
Pad size	46 inches (116.8cm) long x 20 inches (50.8cm) wide	36 inches (91.4cm) long x 20 inches (50.8cm) wide	36 inches (91.4cm) long x 20 inches (50.8cm) wide
Pad thickness	1/2-inch (1.25cm)	1/4-inch (0.63cm)	1/8-inch (0.31cm)
Key features	Pressure reduction	Patient comfort, lightweight, pressure reduction	Versatile and lightweight
Patient size	25lbs. (11.3kg) and up	0.8 lbs. (0.36kg) and up	0.8 lbs. (0.36kg) and up
Pad weight	18.5 lbs. (53.7kg)	8 lbs. (3.62kg)	4 lbs. (1.81kg)
Pad use	One-sided	One-sided	Two-sided
Lightweight		✓	✓
Flexible		✓	✓
Pressure reduction	✓	✓	
Comfort	✓	✓	
Max life	24 months	24 months	24 months
Additional information	Single cord option. Dual cord option for use with two generators.	Single cord option. Dual cord option for use with two generators. Pediatrics to adults.	Single cord option. Dual cord option for use with two generators. Pediatrics to adults.

1. (C2476)

Please refer always to the Instructions for Use / Package Insert that come with the device for the most current and complete instructions.

Ethicon Endo-Surgery (Europe) GmbH
Hummelsbütteler Steindamm 71
22851 Norderstedt, Germany

www.ethicon.com