

Reload comparisons based on extensive analysis of closed staple heights

Endo GIA™ with Tri-Staple™ Technology Per internal testing			ECHELON FLEX™ GST System ¹ Per IFU				
.88-1.5mm	1.5-2.25mm	2.25-3.0mm	1.0-2.0mm	1.5-2.4mm	1.8-3.0mm	2.0-3.3mm	2.3-4.0mm
INTENDED TISSUE THICKNESS Per manufacturer IFU			TISSUE THICKNESS RANGE Porcine tissue measured at 8g/mm ² prior to firing				

¹ System components include ECHELON FLEX™ Powered Plus Stapler and ENDOPATH ECHELON™ Reloads with Gripping Surface Technology.

www.ethicon.com

©2016 Ethicon, Inc. All rights reserved. 060421161213

The third-party trademarks used herein are trademarks of their respective owners.

For complete product details, see Instructions for Use.