


ENDOPATH® Probe Plus II

Multi-Functional Electrosurgery System



Pistol grip handle now requires
**73% less force to activate
suction and irrigation***

Cautery, suction and irrigation in a single device

The ENDOPATH® Probe Plus II Multi-Functional Electrosurgery System offers **two handles** and **seven electrode/shaft options** to meet a wide variety of surgical needs.

Features

- Performs dissection, coagulation, fluid evacuation and irrigation
- Improved pistol grip handle for optimized ergonomics
- 5mm shaft diameter can be used with a variety of trocar sizes
- Compatible with commonly used irrigation and suction devices, as well as monopolar electrosurgery units
- Interchangeable shafts allow for multiple configurations, eliminating the need for additional devices

Talk to your Ethicon representative to learn more.

Top view of pistol grip



Pistol grip handle



Pencil grip handle



Electrode options



Hook



Curved
Dissector



Spatula



Right Angle

Shafts allow monopolar electrosurgery with suction and irrigation

ENDOPATH® Probe Plus II

Multi-Functional Electrosurgery System

Ordering information

Product code	Description	Shaft diameter (mm)	Shaft length (cm)	Instruments per sales unit
Handles				
EPHO2	Hand Controlled, Pistol Grip Handle			6
EPHO4	Hand Controlled, Pencil Grip Handle			6
Electrodes				
EPSO1	Hook Electrode, Hollow Tip	5mm	34cm	6
EPSO2	Spatula Electrode, Hollow Tip	5mm	34cm	6
EPSO3	Right Angle Electrode, Hollow Tip	5mm	34cm	6
EPSO4	Curved Dissector Electrode, Hollow Tip	5mm	34cm	6
EPSO5	Hook Electrode, Hollow Tip	5mm	29cm	6
EPSO6	Spatula Electrode, Hollow Tip	5mm	29cm	6
EPSO7	Right Angle Electrode, Hollow Tip	5mm	29cm	6
Tubing				
EPTO1	Dual-Bottle Irrigation Tube Set			6
EPTO3	Dual-Bottle Irrigation Tube Set			6

How to order

All ENDOPATH® Probe Plus II System purchase orders are made to Johnson & Johnson Healthcare Carte Systems, Inc. (JJHCS).

Electronic ordering options

Note: Placing order electronically avoids minimum order fees for hospitals.

Johnson & Johnson Gateway

Visit jnsgateway.com or call **1-866-JNJ-GATE**.

Global Healthcare Exchange

Visit ghx.com or call **1-800-YOUR-GHX**.

Electronic Data Interchange

Call the JJHCS Help Line at **1-800-262-2888**.

Customer support

For more information or product support, call **1-877-ETHICON** or visit ethicon.com/healthcare-professionals/products/other/electrosurgical-irrigation.

See **Instructions for Use** for complete product details.

Non-electronic/manual ordering options

Call JJHCS at **1-800-255-2500** (option 1) between 8:30 am and 8:00 pm Eastern time or fax your order to **1-732-562-2212**.