

With the Ethicon commitment to improved surgical experiences the difference is in the details.

ENDOPOUCH® Specimen Retrieval Bag

Ethicon believes every device used in the OR plays a critical role. That is why we work to ensure that each device is engineered to exacting standards and delivers performance you expect. Consider the important differences in design and ease of use between the ENDOPOUCH Specimen Retrieval Bag compared to Inzii™ Retrieval System by Applied Medical.



ENDOPOUCH Specimen Retrieval Bag:
requires on average
92% less peak force to deploy
(1.55 lbf compared to 20.21 lbf)

ENDOPOUCH Specimen Retrieval Bag:
27% greater bag length
(6.44 in compared to 5.08 in)



ENDOPOUCH Specimen Retrieval Bag:
12% smaller grip span
(4.00 in compared to 4.56 in)

Ease of use and proven performance*

A smaller, more comfortable grip span gives the ENDOPOUCH Specimen Retrieval Bag an ergonomic advantage for easier one-handed operation. This, combined with a low insertion force and smooth deployment, contributes to a comfortable and intuitive surgical experience.

At Ethicon, we take pride in all of our surgical devices — and their critical impact on human lives. The ENDOPOUCH Specimen Retrieval Bag is evidence of this ongoing commitment.