

# Advanced Bipolar Energy Quick Comparison

	ENSEAL® G2 Tissue Sealers				LigaSure™			Olympus
	Curved Jaw	Straight Jaw	Articulating Curved Jaw	Articulating Straight Jaw	5mm Blunt Tip	Advance	Maryland Jaw	THUNDERBEAT
			 Articulating Curved Jaw shown					
① Shaft Diameter (mm)	5	5	5	5	5	5	5	5
② Cut Length (mm)	17	16	17	16	16	16	18	16
③ Jaw Length (mm)	19	20	20	20	20	18	21	17
④ Jaw Aperture (mm)	17	15	18	15	12	15	14	13
⑤ Jaw Width (mm)	3	5	3	5	5	4	4	3
⑥ Shaft Length (cm)	14, 25, 35, 45	14, 25, 35, 45	35, 45	35, 45	20, 37, 44	44	23, 37, 44	10, 20, 35, 45