

**SELF-DRILLING  
SPEED**


**MatrixRIB®**

FIXATION SYSTEM

# Improved procedural efficiency<sup>1</sup>




10 mm



20 mm

MatrixRIB® Fixation System Self-Drilling Screws improve procedural efficiency by reducing the number of steps required to insert screws compared to self-tapping screws.



## Reliable screw insertion<sup>2</sup>

Screw Guides ensure proper screw insertion angulation.

# SELF-DRILLING SPEED


## Improved locking strength<sup>3</sup>

**3X**

The locking strength of the MatrixRIB Self-Drilling Locking Screw is 3X greater than the Zimmer Biomet self-drilling screws based on cantilever bending tests. \*



Screw Cantilever Bending Load



1. DePuy Synthes Document #0000284579
2. DePuy Synthes Document #0000279377†
3. DePuy Synthes Document #0000280454†

\* Bench test data are not indicative of clinical test results

†DePuy Synthes documents (0000277013 and 0000280454). Bench test, n=6, comparing 20mm MatrixRIB SD locking vs 16mm Zimmer Biomet self-drilling locking screws. Bench testing may not be indicative of clinical performance."

Some devices listed in this brochure may not have been licensed in accordance with Canadian law and may not be for sale in Canada. Please contact your Sales Consultant for items approved for sale in Canada.

Please contact your DePuy Synthes Sales Consultants for details. To order, call DePuy Synthes Customer Service at (800) 523-0322.

Limited Warranty and Disclaimer: DePuy Synthes CMF products are sold with a limited warranty to the original purchaser against defects in workmanship and materials. Any other express or implied warranties, including warranties of merchantability or fitness, are hereby disclaimed.

WARNING: In the USA, this product has labeling limitations. See package insert for complete information.

CAUTION: USA Law restricts these devices to sale by or on the order of a physician.

Not all products are currently available in all markets.

The third party trademarks used herein are the trademarks of their respective owners.



**Synthes USA Products, LLC**  
1302 Wrights Lane East  
West Chester, PA 19380  
Telephone: (610) 719-5000  
To order: (800) 523-0322

**Synthes USA, LLC**  
1101 Synthes Avenue  
Monument, CO 80132

[www.depuysynthes.com](http://www.depuysynthes.com)

©DePuy Synthes 2019. All rights reserved.  
109086-190308 DSUS 05/19

## Adaptive torque limiting Knows when to stop

Patented torque limiting technology detects when a screw is seated and cuts power to the motor to prevent stripping.

### MatrixPRO® II Driver



**90°**

### Minimally invasive 90° intracollers

Multiple attachments for both driving and drilling

- Straight Hex collet
- Straight J-Latch collet
- 90° intracollers, short
- 90° intracollers, long