


# PERMATAPE™ Suture

**Your Technique,  
Your Choice**




**PERMATAPE™ Suture now offers:**


- An optimized solution for both knotless and knot-tying repairs
- A wider contact area with lower abrasiveness and less bulk as compared to the leading competitor<sup>1,2,3</sup>
- Compatibility with most existing instrumentation

# PERMATAPE™ Suture

**100% Ultra-High Molecular Weight Polyethylene (UHMWPE)  
flat braid for strength and handling.**



**76% less abrasive than FiberTape<sup>1</sup>  
30-cycle test with 15N load**



**12% more contact area with  
54% less bulk as compared  
to FiberTape<sup>2,3</sup>**

## Ordering Information

Product	Description
223169	PERMATAPE Suture Blue Braided Flat Suture, 2.5 mm (Single Pack)
223170	PERMATAPE Suture White/Blue Braided Flat Suture, 2.5 mm (Single Pack)
223165	PERMATAPE Suture Blue Braided Flat Suture, 2.5 mm (6 Pack)
223166	PERMATAPE Suture White/Blue Braided Flat Suture, 2.5 mm (6 Pack)
288233	EXPRESSEW III AutoCapture+ Flexible Suture Passer
214141	EXPRESSEW III Needles (5 Pack)
288234	EXPRESSEW III AutoCapture+ Retention Plate with Loader (5 Pack)
288235	EXPRESSEW III AutoCapture+ Retention Plate with Loader (1 Each)



**Compatible with EXPRESSEW® III  
Autocapture+ Suture Passer**

## References

1. 103288856 for data on file
2. 103313374 for data on file
3. 103313375 for data on file

FiberTape is a trademark of Arthrex, Inc.

Not all products are currently available in all markets. The third party trademarks used herein are trademarks of their respective owners.



**Johnson & Johnson Medical Pty Ltd t/a DePuy Synthes**  
1-5 Khartoum Road,  
North Ryde, NSW, 2113  
Australia

**Johnson & Johnson (NZ) Limited**  
507 Mt Wellington Hwy,  
Mt Wellington, Auckland, 1060  
New Zealand