

EVERPOINT® Cardiovascular Needle

Product Catalog

EVERPOINT Needles provide cardiovascular surgeons with a stronger, more bend-resistant needle that stays sharp, pass after pass

Compared to conventional stainless steel needles, EVERPOINT Needles offer^{1,2}:

- **Superior penetration**
- **Up to 121% more bend resistance**
- **Up to 38% more strength**

Patients undergoing CABG surgery often present with comorbidities that create suturing challenges, such as calcified and friable vessels. EVERPOINT Needles are specifically designed to help surgeons avoid unnecessary complications in difficult CABG cases.

Enhanced sharpness

Unique multilayer silicone coating provides **outstanding penetration** and maintains tip durability, even through unexpected calcium deposits^{1,2}

Improved strength and bend resistance


Without compromising ductility, the ultra high-strength **tungsten-rhenium alloy is significantly stronger and more bend resistant** than stainless steel needles, allowing improved handling and placement²

Unsurpassed performance

This unique combination of exceptional strength and durable penetration is what puts EVERPOINT Needles in a class of their own



EVERPOINT® Cardiovascular Needles

a stronger, more bend-resistant needle that stays sharp, pass after pass^{1,2}











EVERPOINT Needle Product Codes







Taper Point ⊙

Needle	Product Code Quantity/Pack x Length			
	8-0	7-0	6-0	5-0
RB-2  18 mil, 13 mm 1/2 Circle			EP8714H 1 x 45 cm / 18"	
			EP8711H 1 x 75 cm / 30"	
			EPM8711 4 x 75 cm / 30"	
RB-2  10.5 mil, 13 mm 1/2 Circle				*EPH8710H 1 x 75 cm / 30"
				*EPH8716H 1 x 90 cm / 36"

*HEMO-SEAL™ Technology

Needle	Product Code Quantity/Pack x Length			
	8-0	7-0	6-0	5-0
 10.5 mil, 13 mm 3/8 Circle		EP8800H 1 x 45 cm / 18"	EP8718H 1 x 45 cm / 18"	*EPH7477H 1 x 75 cm / 30"
		EPM8800 4 x 45 cm / 18"	EPM8718 4 x 45 cm / 18"	*EPH8720H 1 x 90 cm / 36"
		EPM8727 4 x 60 cm / 24"	EP8726H 1 x 60 cm / 24"	*EPH8725H 1 x 60 cm / 24"
		EP8708H 1 x 75 cm / 30"	EPM8726 4 x 60 cm / 24"	*EPH7477H 1 x 75 cm / 30"
			EP8706H 1 x 75 cm / 30"	
			EPM8206 2 x 75 cm / 30"	
 10.5 mil, 13 mm 3/8 Circle Single Armed			EP8889H 1 x 75 cm / 30"	*EPH8890H 1 x 75 cm / 30"
 10 mil, 11 mm 3/8 Circle			EP8610H 1 x 60 cm / 24"	
			EPM8610 4 x 60 cm / 24"	
			EP8776H 1 x 75 cm / 30"	
			EPM8776 4 x 75 cm / 30"	
 10.5 mil, 9.3 mm 3/8 Circle			EP8806H 1 x 45 cm / 18"	
			EP8805H 1 x 60 cm / 24"	
 10 mil, 9.3 mm 3/8 Circle		EPM8203 2 x 75 cm / 30"	EPM8806 4 x 45 cm / 18"	*EPH9702H 1 x 60 cm / 24"
		EPM8703 4 x 75 cm / 30"	EPM8205 2 x 60 cm / 24"	
		EPM8702 4 x 60 cm / 24"	EPM8805 4 x 60 cm / 24"	
			EP8709H 1 x 75 cm / 30"	
			EPM8709 4 x 75 cm / 30"	

Needle	Product Code Quantity/Pack x Length			
	8-0	7-0	6-0	5-0
BV175-8  7.8 mil, 9.3 mm 3/8 Circle	EP8752H 1 x 45 cm / 18"	EP8746H 1 x 45 cm / 18"		
	EPM8752 2 x 45 cm / 18"	EPM8746 2 x 45 cm / 18"		
	EP8753H 1 x 60 cm / 24"	EP8747H 1 x 60 cm / 24"		
	EPM8753 2 x 60 cm / 24"	EPM8747 1 x 60 cm / 24"		
	EP8754H 1 x 75 cm / 30"	EPM8745 4 x 60 cm / 24"		
	EPM8754 2 x 75 cm / 30"	EP8755H 1 x 75 cm / 30"		
		EPM8755 2 x 75 cm / 30"		
		EPM8751 4 x 75 cm / 30"		
BV175-7  7.8 mil, 8 mm 3/8 Circle	EPM8767 2 x 60 cm / 24"	EPM8765 2 x 45 cm / 18"		
		EP8766H 1 x 60 cm / 24"		
		EPM8766 2 x 60 cm / 24"		
		EPM8775 4 x 60 cm / 24"		
		EPM8768 2 x 75 cm / 30"		
		EPM8785 4 x 75 cm / 30"		
BV175-6  8 mil, 8 mm 3/8 Circle	EP8740H 1 x 45 cm / 18"	EP8734H 1 x 45 cm / 18"		
	EPM8742 2 x 45 cm / 18"	EPM8736 2 x 45 cm / 18"		
	EP8741H 1 x 60 cm / 24"	EPM8734 4 x 45 cm / 18"		
	EPM8743 2 x 60 cm / 24"	EP8735H 1 x 60 cm / 24"		
	EPM8741 4 x 60 cm / 24"	EPM8737 2 x 60 cm / 24"		
		EPM8735 4 x 60 cm / 24"		
		EPM8733 2 x 75 cm / 30"		
		EPM8731 4 x 75 cm / 30"		
BV130-5  6 mil, 6.5 mm 3/8 Circle	EP8730H 1 x 45 cm / 18"			
	EPM8738 2 x 45 cm / 18"			
	EPM8730 4 x 45 cm / 18"			
	EP8732H 1 x 60 cm / 24"			
	EPM8739 2 x 60 cm / 24"			
	EPM8732 4 x 60 cm / 24"			

Needle	Product Code Quantity/Pack x Length			
	8-0	7-0	6-0	5-0
CC  11 mil, 9.3 mm 3/8 Circle		EPM8704 4 x 60 cm / 24"		
CC-1  12 mil, 13 mm 3/8 Circle			EP8719H 1 x 45 cm / 18" EPM8719 4 x 45 cm / 18" EPM8707 4 x 75 cm / 30"	
CC-1  10 mil, 13 mm 3/8 Circle			EP8707H 1 x 75 cm / 30"	
CC  10 mil, 9.3 mm 3/8 Circle		EPM8705 4 x 45 cm / 18"	EP8807H 1 x 60 cm / 24"	
CC  9.1 mil, 9.3 mm 3/8 Circle			EPM8807 4 x 60 cm / 24"	
CC175-8  8 mil, 9.3 mm 3/8 Circle		EPM8704SL 4 x 60 cm / 24" EPM8812SL 4 x 75 cm / 30"		

To order EVERPOINT® Cardiovascular Needles today, contact your Ethicon Representative or call 1-877-ETHICON (384-4266).

For complete indications, contraindications, warnings, precautions, and adverse reactions, please reference full package insert.

References

1. EVERPOINT Competitive Assessment, CT11-013, 2011. Ethicon, Inc. 2. Cichocki FR Jr, Maurer RE, Bar SN. Tungsten-rhenium suture needles with improved properties for coronary artery bypass graft surgery. *J Biomed Mater Res B Appl Biomater*. 2010;94(2):493-500.