

# INVENTORY CONTROL FORM



## Universal Small Fragment System

PART OF THE *Johnson & Johnson* FAMILY OF COMPANIES

3.5 mm VA LCP® Proximal Tibia Implant Set 01.133.213

**Case ID:**

**To Order (USA): (800) 523-0322**  
**To Order (Canada): (844) 243-4321**

Date: \_\_\_\_\_

Hospital: \_\_\_\_\_

Surgeon: \_\_\_\_\_

Procedure: \_\_\_\_\_

## Implants

### 3.5 mm VA LCP® Proximal Tibia Plate, Small Bend

Available nonsterile or sterile-packed. Add "S" to catalog number to order sterile product.

<input type="checkbox"/>	<b>02.127.210</b>	4 holes, 87 mm, right
<input type="checkbox"/>	<b>02.127.211</b>	4 holes, 87 mm, left
<input type="checkbox"/>	<b>02.127.220</b>	6 holes, 117 mm, right
<input type="checkbox"/>	<b>02.127.221</b>	6 holes, 117 mm, left
<input type="checkbox"/>	<b>02.127.230</b>	8 holes, 147 mm, right
<input type="checkbox"/>	<b>02.127.231</b>	8 holes, 147 mm, left
<input type="checkbox"/>	<b>02.127.240</b>	10 holes, 177 mm, right
<input type="checkbox"/>	<b>02.127.241</b>	10 holes, 177 mm, left



### 3.5 mm VA LCP® Proximal Tibia Plate, Large Bend

Available nonsterile or sterile-packed. Add "S" to catalog number to order sterile product.

<input type="checkbox"/>	<b>02.127.310</b>	4 holes, 87 mm, right
<input type="checkbox"/>	<b>02.127.311</b>	4 holes, 87 mm, left
<input type="checkbox"/>	<b>02.127.320</b>	6 holes, 117 mm, right
<input type="checkbox"/>	<b>02.127.321</b>	6 holes, 117 mm, left
<input type="checkbox"/>	<b>02.127.330</b>	8 holes, 147 mm, right
<input type="checkbox"/>	<b>02.127.331</b>	8 holes, 147 mm, left
<input type="checkbox"/>	<b>02.127.340</b>	10 holes, 177 mm, right
<input type="checkbox"/>	<b>02.127.341</b>	10 holes, 177 mm, left



### 3.5 mm Variable Angle Locking Screw, Self-Tapping, T15 StarDrive™ Recess



<input type="checkbox"/>	<b>02.127.110</b>	10 mm	<input type="checkbox"/>	<b>02.127.136</b>	36 mm
<input type="checkbox"/>	<b>02.127.112</b>	12 mm	<input type="checkbox"/>	<b>02.127.138</b>	38 mm
<input type="checkbox"/>	<b>02.127.114</b>	14 mm	<input type="checkbox"/>	<b>02.127.140</b>	40 mm
<input type="checkbox"/>	<b>02.127.116</b>	16 mm	<input type="checkbox"/>	<b>02.127.142</b>	42 mm
<input type="checkbox"/>	<b>02.127.118</b>	18 mm	<input type="checkbox"/>	<b>02.127.144</b>	44 mm
<input type="checkbox"/>	<b>02.127.120</b>	20 mm	<input type="checkbox"/>	<b>02.127.146</b>	46 mm
<input type="checkbox"/>	<b>02.127.122</b>	22 mm	<input type="checkbox"/>	<b>02.127.148</b>	48 mm
<input type="checkbox"/>	<b>02.127.124</b>	24 mm	<input type="checkbox"/>	<b>02.127.150</b>	50 mm
<input type="checkbox"/>	<b>02.127.126</b>	26 mm	<input type="checkbox"/>	<b>02.127.152</b>	52 mm
<input type="checkbox"/>	<b>02.127.128</b>	28 mm	<input type="checkbox"/>	<b>02.127.154</b>	54 mm
<input type="checkbox"/>	<b>02.127.130</b>	30 mm	<input type="checkbox"/>	<b>02.127.156</b>	56 mm
<input type="checkbox"/>	<b>02.127.132</b>	32 mm	<input type="checkbox"/>	<b>02.127.158</b>	58 mm
<input type="checkbox"/>	<b>02.127.134</b>	34 mm	<input type="checkbox"/>	<b>02.127.160</b>	60 mm

The following screw sizes may be stored in the screw rack of the VA LCP Proximal Tibia Implant Set

<input type="checkbox"/>	<b>02.127.165</b>	65 mm*	<input type="checkbox"/>	<b>02.127.180</b>	80 mm*
<input type="checkbox"/>	<b>02.127.170</b>	70 mm*	<input type="checkbox"/>	<b>02.127.185</b>	85 mm*
<input type="checkbox"/>	<b>02.127.175</b>	75 mm*	<input type="checkbox"/>	<b>02.127.190</b>	90 mm*

\* 2.7 mm screws longer than 68 mm and 3.5 mm screws longer than 60 mm will not fit in the screw rack of the Universal Small Fragment System.

# INVENTORY CONTROL FORM

## Universal Small Fragment System

3.5 mm VA LCP® Proximal Tibia Implant Set 01.133.213

### Implants (continued)

#### 3.5 mm Locking Screw, Self-Tapping, T15 StarDrive™ Recess



The following screw sizes may be stored in the screw rack of the Universal Small Fragment System

<input type="checkbox"/>	<b>212.101</b>	10 mm	<input type="checkbox"/>	<b>212.116</b>	38 mm
<input type="checkbox"/>	<b>212.102</b>	12 mm	<input type="checkbox"/>	<b>212.117</b>	40 mm
<input type="checkbox"/>	<b>212.103</b>	14 mm	<input type="checkbox"/>	<b>212.118</b>	42 mm
<input type="checkbox"/>	<b>212.104</b>	16 mm	<input type="checkbox"/>	<b>212.134</b>	44 mm
<input type="checkbox"/>	<b>212.105</b>	18 mm	<input type="checkbox"/>	<b>212.119</b>	45 mm
<input type="checkbox"/>	<b>212.106</b>	20 mm	<input type="checkbox"/>	<b>212.136</b>	46 mm
<input type="checkbox"/>	<b>212.107</b>	22 mm	<input type="checkbox"/>	<b>212.120</b>	48 mm
<input type="checkbox"/>	<b>212.108</b>	24 mm	<input type="checkbox"/>	<b>212.121</b>	50 mm
<input type="checkbox"/>	<b>212.109</b>	26 mm	<input type="checkbox"/>	<b>212.122</b>	52 mm
<input type="checkbox"/>	<b>212.110</b>	28 mm	<input type="checkbox"/>	<b>02.212.054</b>	54 mm
<input type="checkbox"/>	<b>212.111</b>	30 mm	<input type="checkbox"/>	<b>212.123</b>	55 mm
<input type="checkbox"/>	<b>212.112</b>	32 mm	<input type="checkbox"/>	<b>02.212.056</b>	56 mm
<input type="checkbox"/>	<b>212.113</b>	34 mm	<input type="checkbox"/>	<b>02.212.058</b>	58 mm
<input type="checkbox"/>	<b>212.115</b>	36 mm	<input type="checkbox"/>	<b>212.124</b>	60 mm

The following screw sizes may be stored in the screw rack of the VA LCP Proximal Tibia Implant Set

<input type="checkbox"/>	<b>212.125</b>	65 mm*	<input type="checkbox"/>	<b>212.128</b>	80 mm*
<input type="checkbox"/>	<b>212.126</b>	70 mm*	<input type="checkbox"/>	<b>212.129</b>	85 mm*
<input type="checkbox"/>	<b>212.127</b>	75 mm*	<input type="checkbox"/>	<b>212.130</b>	90 mm*

#### 3.5 mm Cortex Screw, 2.5 mm Hex Recess



The following screw sizes may be stored in the screw rack of the Universal Small Fragment System

<input type="checkbox"/>	<b>204.810</b>	10 mm	<input type="checkbox"/>	<b>204.834</b>	34 mm
<input type="checkbox"/>	<b>204.812</b>	12 mm	<input type="checkbox"/>	<b>204.836</b>	36 mm
<input type="checkbox"/>	<b>204.814</b>	14 mm	<input type="checkbox"/>	<b>204.838</b>	38 mm
<input type="checkbox"/>	<b>204.816</b>	16 mm	<input type="checkbox"/>	<b>204.840</b>	40 mm
<input type="checkbox"/>	<b>204.818</b>	18 mm	<input type="checkbox"/>	<b>204.842</b>	42 mm
<input type="checkbox"/>	<b>204.820</b>	20 mm	<input type="checkbox"/>	<b>204.844</b>	44 mm
<input type="checkbox"/>	<b>204.822</b>	22 mm	<input type="checkbox"/>	<b>204.845</b>	45 mm
<input type="checkbox"/>	<b>204.824</b>	24 mm	<input type="checkbox"/>	<b>204.846</b>	46 mm
<input type="checkbox"/>	<b>204.826</b>	26 mm	<input type="checkbox"/>	<b>204.848</b>	48 mm
<input type="checkbox"/>	<b>204.828</b>	28 mm	<input type="checkbox"/>	<b>204.850</b>	50 mm
<input type="checkbox"/>	<b>204.830</b>	30 mm	<input type="checkbox"/>	<b>204.855</b>	55 mm
<input type="checkbox"/>	<b>204.832</b>	32 mm	<input type="checkbox"/>	<b>204.860</b>	60 mm






The following screw sizes may be stored in the screw rack of the VA LCP Proximal Tibia Implant Set

<input type="checkbox"/>	<b>204.865</b>	65 mm*	<input type="checkbox"/>	<b>204.880</b>	80 mm*
<input type="checkbox"/>	<b>204.870</b>	70 mm*	<input type="checkbox"/>	<b>204.885</b>	85 mm*
<input type="checkbox"/>	<b>204.875</b>	75 mm*	<input type="checkbox"/>	<b>204.890</b>	90 mm*

### Instruments and Other Items




#### Drill Bits

The following drill bits are contained in the Universal Small Fragment Insertion Tray  
Available nonsterile or sterile-packed. Add "S" to catalog number to order sterile product.


<input type="checkbox"/>	<b>03.133.102</b>	2.5 mm Drill Bit, Quick Coupling, 45 mm	
<input type="checkbox"/>	<b>03.133.103</b>	2.5 mm Drill Bit, Quick Coupling, 80 mm	
<input type="checkbox"/>	<b>03.133.106</b>	2.8 mm Drill Bit, Quick Coupling, 45 mm	
<input type="checkbox"/>	<b>03.133.107</b>	2.8 mm Drill Bit, Quick Coupling, 80 mm	
<input type="checkbox"/>	<b>03.133.109</b>	3.5 mm Drill Bit, Quick Coupling, 60 mm	

The following drill bits are contained in the Universal Small Fragment VA LCP Proximal Tibia Implant Anatomy Tray

Available nonsterile or sterile-packed. Add "S" to catalog number to order sterile product.

<input type="checkbox"/>	<b>03.133.104</b>	2.5 mm Drill Bit, Quick Coupling, 150 mm	
<input type="checkbox"/>	<b>03.133.108</b>	2.8 mm Drill Bit, Quick Coupling, 110 mm	
<input type="checkbox"/>	<b>03.133.110</b>	3.5 mm Drill Bit, Quick Coupling, 105 mm	

#### Other Items

<input type="checkbox"/>	<b>292.20</b>	2.0 mm Kirschner Wire with trocar point, 150 mm	
--------------------------	---------------	---	---

\*2.7 mm screws longer than 68 mm and 3.5 mm screws longer than 60 mm will not fit in the screw rack of the Universal Small Fragment System.