

# Inventory Control Form

## 4.5, 5.5, 6.5, 7.5 mm Cannulated Compression Headless Screws (CCHS)

Case Information (eg, Case ID):

DO NOT PUT ANY  
PATIENT IDENTIFIABLE  
INFORMATION IN  
THIS FIELD

To order (USA): 800-523-0322

To order (Canada): 844-243-4321

Date: \_\_\_\_\_

Hospital: \_\_\_\_\_

Surgeon: \_\_\_\_\_

Procedure: \_\_\_\_\_

### Implants The screws are made of titanium alloy (TAN)

#### 4.5 mm Cannulated Compression Headless Screw



Length	Short Thread Part Number	Length	Long Thread Part Number
		N/A	
		N/A	
		N/A	
		N/A	
		N/A	
		30 mm	04.334.530
		32 mm	04.334.532
		34 mm	04.334.534
		36 mm	04.334.536
		38 mm	04.334.538
		40 mm	04.334.540
		42 mm	04.334.542
		44 mm	04.334.544
		46 mm	04.334.546
		48 mm	04.334.548
		50 mm	04.334.550
		55 mm	04.334.555
		60 mm	04.334.560
		65 mm	04.334.565
		70 mm	04.334.570
		75 mm	04.334.575
		80 mm	04.334.580
		85 mm	04.334.585
		90 mm	04.334.590
		95 mm	04.334.595
		100 mm	04.334.500
		105 mm	04.334.501
		110 mm	04.334.502

Thread Length  
Short Thread: approx. 25% Total Length or 8 mm  
Long Thread: approx. 40% Total Length

#### 5.5 mm Cannulated Compression Headless Screw



Length	Short Thread Part Number	Length	Long Thread Part Number
		N/A	
		N/A	
		N/A	
		N/A	
		N/A	
		30 mm	04.334.630
		32 mm	04.334.632
		34 mm	04.334.634
		36 mm	04.334.636
		38 mm	04.334.638
		40 mm	04.334.640
		42 mm	04.334.642
		44 mm	04.334.644
		46 mm	04.334.646
		48 mm	04.334.648
		50 mm	04.334.650
		55 mm	04.334.655
		60 mm	04.334.660
		65 mm	04.334.665
		70 mm	04.334.670
		75 mm	04.334.675
		80 mm	04.334.680
		85 mm	04.334.685
		90 mm	04.334.690
		95 mm	04.334.695
		100 mm	04.334.600
		105 mm	04.334.601
		110 mm	04.334.602

Thread Length  
Short Thread: approx. 25% Total Length or 8 mm  
Long Thread: approx. 40% Total Length

Available packed in sterile tubes. Add suffix "TS" to part number. Sterile part availability might be different depending on the country or region.

# Inventory Control Form

## 4.5, 5.5, 6.5, 7.5 mm Cannulated Compression Headless Screws (CCHS)

### Implants The screws are made of titanium alloy (TAN)

#### 6.5 mm Cannulated Compression Headless Screw <sup>◇</sup>



Length	Short Thread Part Number	Length	Long Thread Part Number
<input type="checkbox"/> 30 mm	04.333.730	N/A	
<input type="checkbox"/> 35 mm	04.333.735	N/A	
<input type="checkbox"/> 40 mm	04.333.740	N/A	
<input type="checkbox"/> 45 mm	04.333.745	<input type="checkbox"/> 45 mm	04.334.745
<input type="checkbox"/> 50 mm	04.333.750	<input type="checkbox"/> 50 mm	04.334.750
<input type="checkbox"/> 55 mm	04.333.755	<input type="checkbox"/> 55 mm	04.334.755
<input type="checkbox"/> 60 mm	04.333.760	<input type="checkbox"/> 60 mm	04.334.760
<input type="checkbox"/> 65 mm	04.333.765	<input type="checkbox"/> 65 mm	04.334.765
<input type="checkbox"/> 70 mm	04.333.770	<input type="checkbox"/> 70 mm	04.334.770
<input type="checkbox"/> 75 mm	04.333.775	<input type="checkbox"/> 75 mm	04.334.775
<input type="checkbox"/> 80 mm	04.333.780	<input type="checkbox"/> 80 mm	04.334.780
<input type="checkbox"/> 85 mm	04.333.785	<input type="checkbox"/> 85 mm	04.334.785
<input type="checkbox"/> 90 mm	04.333.790	<input type="checkbox"/> 90 mm	04.334.790
<input type="checkbox"/> 95 mm	04.333.795	<input type="checkbox"/> 95 mm	04.334.795
<input type="checkbox"/> 100 mm	04.333.700	<input type="checkbox"/> 100 mm	04.334.700
<input type="checkbox"/> 105 mm	04.333.701	<input type="checkbox"/> 105 mm	04.334.701
<input type="checkbox"/> 110 mm	04.333.702	<input type="checkbox"/> 110 mm	04.334.702
<input type="checkbox"/> 115 mm	04.333.703	<input type="checkbox"/> 115 mm	04.334.703
<input type="checkbox"/> 120 mm	04.333.704	<input type="checkbox"/> 120 mm	04.334.704
<input type="checkbox"/> 125 mm	04.333.705	<input type="checkbox"/> 125 mm	04.334.705
<input type="checkbox"/> 130 mm	04.333.706	<input type="checkbox"/> 130 mm	04.334.706

#### 7.5 mm Cannulated Compression Headless Screw <sup>◇</sup>



Length	Short Thread Part Number	Length	Long Thread Part Number
<input type="checkbox"/> 30 mm	04.333.830	N/A	
<input type="checkbox"/> 35 mm	04.333.835	N/A	
<input type="checkbox"/> 40 mm	04.333.840	N/A	
<input type="checkbox"/> 45 mm	04.333.845	<input type="checkbox"/> 45 mm	04.334.845
<input type="checkbox"/> 50 mm	04.333.850	<input type="checkbox"/> 50 mm	04.334.850
<input type="checkbox"/> 55 mm	04.333.855	<input type="checkbox"/> 55 mm	04.334.855
<input type="checkbox"/> 60 mm	04.333.860	<input type="checkbox"/> 60 mm	04.334.860
<input type="checkbox"/> 65 mm	04.333.865	<input type="checkbox"/> 65 mm	04.334.865
<input type="checkbox"/> 70 mm	04.333.870	<input type="checkbox"/> 70 mm	04.334.870
<input type="checkbox"/> 75 mm	04.333.875	<input type="checkbox"/> 75 mm	04.334.875
<input type="checkbox"/> 80 mm	04.333.880	<input type="checkbox"/> 80 mm	04.334.880
<input type="checkbox"/> 85 mm	04.333.885	<input type="checkbox"/> 85 mm	04.334.885
<input type="checkbox"/> 90 mm	04.333.890	<input type="checkbox"/> 90 mm	04.334.890
<input type="checkbox"/> 95 mm	04.333.895	<input type="checkbox"/> 95 mm	04.334.895
<input type="checkbox"/> 100 mm	04.333.800	<input type="checkbox"/> 100 mm	04.334.800
<input type="checkbox"/> 105 mm	04.333.801	<input type="checkbox"/> 105 mm	04.334.801
<input type="checkbox"/> 110 mm	04.333.802	<input type="checkbox"/> 110 mm	04.334.802
<input type="checkbox"/> 115 mm	04.333.803	<input type="checkbox"/> 115 mm	04.334.803
<input type="checkbox"/> 120 mm	04.333.804	<input type="checkbox"/> 120 mm	04.334.804
<input type="checkbox"/> 125 mm	04.333.805	<input type="checkbox"/> 125 mm	04.334.805
<input type="checkbox"/> 130 mm	04.333.806	<input type="checkbox"/> 130 mm	04.334.806
<input type="checkbox"/> 135 mm	04.333.807	<input type="checkbox"/> 135 mm	04.334.807
<input type="checkbox"/> 140 mm	04.333.808	<input type="checkbox"/> 140 mm	04.334.808

<sup>◇</sup>Available packed in sterile tubes. Add suffix "TS" to part number.  
Sterile part availability might be different depending on the country or region.

# Inventory Control Form

## 4.5, 5.5, 6.5, 7.5 mm Cannulated Compression Headless Screws (CCHS)

### Instruments

#### Instruments for 4.5/5.5 mm Screws

03.333.003 1.6 mm Guide Wire, L 220 mm, Trocar Tip  
 03.333.003S\* (6 pk)



03.333.103 3.0 mm Cannulated Drill Bit, QC  
 03.333.103S\*



03.333.203 4.5/5.5 mm Cannulated Countersink, QC  
 03.333.203S\*



#### Instruments for 6.5/7.5 mm Screws

03.333.004 2.8 mm Guide Wire, L 220 mm, Fluted Tip  
 03.333.004S\* (6 pk)



03.333.104 5.0 mm Cannulated Drill Bit, QC  
 03.333.104S\*



03.333.204 6.5/7.5 mm Cannulated Countersink, QC  
 03.333.204S\*



\*Part numbers with a trailing 'S' designate the sterile part number.  
 Sterile part availability might be different depending on the country or region.



PART OF THE *Johnson & Johnson* FAMILY OF COMPANIES

Manufactured by:  
**Synthes USA, LLC**  
 1101 Synthes Avenue  
 Monument, CO 80132

Manufactured by:  
**Tyber Medical, LLC**  
 83 South Commerce Way, Suite 310  
 Bethlehem, PA 18017

To order (USA): 800-523-0322  
 To order (Canada): 844-243-4321

[www.tybermedical.com](http://www.tybermedical.com)

[jnjmedicaldevices.com](http://jnjmedicaldevices.com)

