

# Inventory Control Form

## QIS and 2.0, 2.5, 3.0, 3.5, 4.0 mm Cannulated Compression Headless Screws (CCHS)

Case Information (eg, Case ID):

DO NOT PUT ANY  
PATIENT IDENTIFIABLE  
INFORMATION IN  
THIS FIELD

To order (USA): 800-523-0322

To order (Canada): 844-243-4321

Date: \_\_\_\_\_

Hospital: \_\_\_\_\_

Surgeon: \_\_\_\_\_

Procedure: \_\_\_\_\_

### Implants The screws are made of titanium alloy (TAN)

#### 2.0 mm Quick Insertion Screw



Length (mm) Part Number

<input type="checkbox"/>	11 mm	04.228.511
<input type="checkbox"/>	12 mm	04.228.512
<input type="checkbox"/>	13 mm	04.228.513
<input type="checkbox"/>	14 mm	04.228.514
<input type="checkbox"/>	15 mm	04.228.515
<input type="checkbox"/>	16 mm	04.228.516
<input type="checkbox"/>	17 mm	04.228.517
<input type="checkbox"/>	18 mm	04.228.518

#### 2.0 mm Cannulated Compression Headless Screw



Length Short Thread Part Number Length Long Thread Part Number

<input type="checkbox"/>	10 mm	04.333.010	N/A		
<input type="checkbox"/>	12 mm	04.333.012	N/A		
<input type="checkbox"/>	14 mm	04.333.014	N/A		
<input type="checkbox"/>	16 mm	04.333.016	N/A		
<input type="checkbox"/>	18 mm	04.333.018	N/A		
<input type="checkbox"/>	20 mm	04.333.020		20 mm	04.334.020
<input type="checkbox"/>	22 mm	04.333.022		22 mm	04.334.022
<input type="checkbox"/>	24 mm	04.333.024		24 mm	04.334.024
<input type="checkbox"/>	26 mm	04.333.026		26 mm	04.334.026
<input type="checkbox"/>	28 mm	04.333.028		28 mm	04.334.028
<input type="checkbox"/>	30 mm	04.333.030		30 mm	04.334.030

Thread Length  
Short Thread: approx. 25% Total Length or 4 mm  
Long Thread: approx. 40% Total Length

#### 2.5 mm Cannulated Compression Headless Screw



Length Short Thread Part Number Length Long Thread Part Number

<input type="checkbox"/>	10 mm	04.333.110	N/A		
<input type="checkbox"/>	12 mm	04.333.112	N/A		
<input type="checkbox"/>	14 mm	04.333.114	N/A		
<input type="checkbox"/>	16 mm	04.333.116	N/A		
<input type="checkbox"/>	18 mm	04.333.118	N/A		
<input type="checkbox"/>	20 mm	04.333.120		20 mm	04.334.120
<input type="checkbox"/>	22 mm	04.333.122		22 mm	04.334.122
<input type="checkbox"/>	24 mm	04.333.124		24 mm	04.334.124
<input type="checkbox"/>	26 mm	04.333.126		26 mm	04.334.126
<input type="checkbox"/>	28 mm	04.333.128		28 mm	04.334.128
<input type="checkbox"/>	30 mm	04.333.130		30 mm	04.334.130
<input type="checkbox"/>	32 mm	04.333.132		32 mm	04.334.132
<input type="checkbox"/>	34 mm	04.333.134		34 mm	04.334.134
<input type="checkbox"/>	36 mm	04.333.136		36 mm	04.334.136
<input type="checkbox"/>	38 mm	04.333.138		38 mm	04.334.138
<input type="checkbox"/>	40 mm	04.333.140		40 mm	04.334.140

Thread Length  
Short Thread: approx. 25% Total Length or 4 mm  
Long Thread: approx. 40% Total Length

∅All implants are also available packed in sterile tubes. Add suffix "TS" to part number.  
Sterile part availability might be different depending on the country or region.

# Inventory Control Form

## QIS and 2.0, 2.5, 3.0, 3.5, 4.0 mm Cannulated Compression Headless Screws (CCHS)

### Implants The screws are made of titanium alloy (TAN)

#### 3.0 mm Cannulated Compression Headless Screw



Length	Short Thread Part Number	Length	Long Thread Part Number
<input type="checkbox"/>	10 mm	04.333.210	N/A
<input type="checkbox"/>	12 mm	04.333.212	N/A
<input type="checkbox"/>	14 mm	04.333.214	N/A
<input type="checkbox"/>	16 mm	04.333.216	N/A
<input type="checkbox"/>	18 mm	04.333.218	N/A
<input type="checkbox"/>	20 mm	04.333.220	20 mm 04.334.220
<input type="checkbox"/>	22 mm	04.333.222	22 mm 04.334.222
<input type="checkbox"/>	24 mm	04.333.224	24 mm 04.334.224
<input type="checkbox"/>	26 mm	04.333.226	26 mm 04.334.226
<input type="checkbox"/>	28 mm	04.333.228	28 mm 04.334.228
<input type="checkbox"/>	30 mm	04.333.230	30 mm 04.334.230
<input type="checkbox"/>	32 mm	04.333.232	32 mm 04.334.232
<input type="checkbox"/>	34 mm	04.333.234	34 mm 04.334.234
<input type="checkbox"/>	36 mm	04.333.236	36 mm 04.334.236
<input type="checkbox"/>	38 mm	04.333.238	38 mm 04.334.238
<input type="checkbox"/>	40 mm	04.333.240	40 mm 04.334.240

Thread Length  
 Short Thread: approx. 25% Total Length or 4 mm  
 Long Thread: approx. 40% Total Length

#### 3.5 mm Cannulated Compression Headless Screw



Length	Short Thread Part Number	Length	Long Thread Part Number
<input type="checkbox"/>	14 mm	04.333.314	N/A
<input type="checkbox"/>	16 mm	04.333.316	N/A
<input type="checkbox"/>	18 mm	04.333.318	N/A
<input type="checkbox"/>	20 mm	04.333.320	N/A
<input type="checkbox"/>	22 mm	04.333.322	N/A
<input type="checkbox"/>	24 mm	04.333.324	24 mm 04.334.324
<input type="checkbox"/>	26 mm	04.333.326	26 mm 04.334.326
<input type="checkbox"/>	28 mm	04.333.328	28 mm 04.334.328
<input type="checkbox"/>	30 mm	04.333.330	30 mm 04.334.330
<input type="checkbox"/>	32 mm	04.333.332	32 mm 04.334.332
<input type="checkbox"/>	34 mm	04.333.334	34 mm 04.334.334
<input type="checkbox"/>	36 mm	04.333.336	36 mm 04.334.336
<input type="checkbox"/>	38 mm	04.333.338	38 mm 04.334.338
<input type="checkbox"/>	40 mm	04.333.340	40 mm 04.334.340
<input type="checkbox"/>	42 mm	04.333.342	42 mm 04.334.342
<input type="checkbox"/>	44 mm	04.333.344	44 mm 04.334.344
<input type="checkbox"/>	46 mm	04.333.346	46 mm 04.334.346
<input type="checkbox"/>	48 mm	04.333.348	48 mm 04.334.348
<input type="checkbox"/>	50 mm	04.333.350	50 mm 04.334.350

Thread Length  
 Short Thread: approx. 25% Total Length or 4 mm  
 Long Thread: approx. 40% Total Length

◇All implants are also available packed in sterile tubes. Add suffix "TS" to part number.  
 Sterile part availability might be different depending on the country or region.

# Inventory Control Form

## QIS and 2.0, 2.5, 3.0, 3.5, 4.0 mm Cannulated Compression Headless Screws (CCHS)

### Implants The screws are made of titanium alloy (TAN)

#### 4.0 mm Cannulated Compression Headless Screw



Length	Short Thread Part Number	Length	Long Thread Part Number		
<input type="checkbox"/>	14 mm	04.333.414	N/A		
<input type="checkbox"/>	16 mm	04.333.416	N/A		
<input type="checkbox"/>	18 mm	04.333.418	N/A		
<input type="checkbox"/>	20 mm	04.333.420	N/A		
<input type="checkbox"/>	22 mm	04.333.422	N/A		
<input type="checkbox"/>	24 mm	04.333.424	<input type="checkbox"/>	24 mm	04.334.424
<input type="checkbox"/>	26 mm	04.333.426	<input type="checkbox"/>	26 mm	04.334.426
<input type="checkbox"/>	28 mm	04.333.428	<input type="checkbox"/>	28 mm	04.334.428
<input type="checkbox"/>	30 mm	04.333.430	<input type="checkbox"/>	30 mm	04.334.430
<input type="checkbox"/>	32 mm	04.333.432	<input type="checkbox"/>	32 mm	04.334.432
<input type="checkbox"/>	34 mm	04.333.434	<input type="checkbox"/>	34 mm	04.334.434
<input type="checkbox"/>	36 mm	04.333.436	<input type="checkbox"/>	36 mm	04.334.436
<input type="checkbox"/>	38 mm	04.333.438	<input type="checkbox"/>	38 mm	04.334.438
<input type="checkbox"/>	40 mm	04.333.440	<input type="checkbox"/>	40 mm	04.334.440
<input type="checkbox"/>	42 mm	04.333.442	<input type="checkbox"/>	42 mm	04.334.442
<input type="checkbox"/>	44 mm	04.333.444	<input type="checkbox"/>	44 mm	04.334.444
<input type="checkbox"/>	46 mm	04.333.446	<input type="checkbox"/>	46 mm	04.334.446
<input type="checkbox"/>	48 mm	04.333.448	<input type="checkbox"/>	48 mm	04.334.448
<input type="checkbox"/>	50 mm	04.333.450	<input type="checkbox"/>	50 mm	04.334.450
<input type="checkbox"/>	52 mm	04.333.452	<input type="checkbox"/>	52 mm	04.334.452
<input type="checkbox"/>	54 mm	04.333.454	<input type="checkbox"/>	54 mm	04.334.454
<input type="checkbox"/>	56 mm	04.333.456	<input type="checkbox"/>	56 mm	04.334.456
<input type="checkbox"/>	58 mm	04.333.458	<input type="checkbox"/>	58 mm	04.334.458
<input type="checkbox"/>	60 mm	04.333.460	<input type="checkbox"/>	60 mm	04.334.460

Thread Length

Short Thread: approx. 25% Total Length or 5 mm

Long Thread: approx. 40% Total Length

∅All implants are also available packed in sterile tubes. Add suffix "TS" to part number. Sterile part availability might be different depending on the country or region.

### Instruments

#### Instruments for 2.0 mm Screws



03.333.000 0.8 mm Guide Wire, L 100 mm Trocar Tip  
 03.333.000S\* (6 pk)



03.333.100 1.6 mm Cannulated Drill Bit, QC  
 03.333.100S\*



03.333.200 2.0 mm Cannulated Countersink, QC  
 03.333.200S\*



#### Instruments for 2.5/3.0 mm Screws

03.333.001 1.1 mm Guide Wire, L 150 mm, Trocar Tip  
 03.333.001S\* (6 pk)



03.333.101 2.0 mm Cannulated Drill Bit, QC  
 03.333.101S\*



03.333.201 2.5/3.0 mm Cannulated Countersink, QC  
 03.333.201S\*

#### Instruments for 3.5/4.0 mm Screws



03.333.002 1.4 mm Guide Wire, L 150 mm, Trocar Tip  
 03.333.002S\* (6 pk)



03.333.102 2.7 mm Cannulated Drill Bit, QC  
 03.333.102S\*



03.333.202 3.5/4.0 mm Cannulated Countersink, QC  
 03.333.202S\*

\*Part numbers with a trailing 'S' designate the sterile part number. Sterile part availability might be different depending on the country or region.



*Manufactured by:*  
**Synthes USA, LLC**  
1101 Synthes Avenue  
Monument, CO 80132

To order (USA): 800-523-0322  
To order (Canada): 844-243-4321

[ijnmedicaldevices.com](http://ijnmedicaldevices.com)

*Manufactured by:*  
**Tyber Medical, LLC**  
83 South Commerce Way, Suite 310  
Bethlehem, PA 18017

[www.tybermedical.com](http://www.tybermedical.com)