

**Features**

- Universal design for left or right femur

**Material (Nail)**

- Titanium-6% aluminum-4% vanadium alloy

**Material (Locking Polymer)**

- Ultrahigh Molecular Weight Polyethylene

**Curvature**


- 1m Radius of curvature

**Cross-Section**


- Ø9 & Ø10 are round
- Ø11 – Ø14 are fluted

**Distal End Diameter**

- Ø9 – Ø12 Nails (5° & 10°): 12mm with 11.2mm flats
- Ø14 Nails: 14mm with no flats



Nail Sizes	Ø 9, 10, Ø 11, 12	Ø 14
	5° & 10°	5° only
160 mm	Available	Not Available
200 mm	Available	Not Available
240 mm	Available	Not Available
280 mm	Available	Available
300 mm	Available	Available
320 mm	Available	Available
340 mm	Available	Available
360 mm	Available	Available
380 mm	Available	Available
400 mm	Available	Available
420 mm	Available	Available
440 mm	Available	Available
460 mm	Available	Available
480 mm	Available	Available



Distal Bend:  
5° & 10° @ 52mm

