

MatrixPRO[®] II Driver

Adaptive Torque Limiting
Knows when to stop

VERSATILE. INTELLIGENT. ADAPTABLE.

Versatile

Multiple Attachments for both driving and drilling.

- Interchangeable collets are compatible with existing DePuy Synthes hex and j-latch driver blades and drill bits to reduce equipment needed in the OR.



Intelligent

Prevents screw stripping and eliminates the need for hand tightening.

- Patented Adaptive Torque Limiting technology detects when a screw is seated and cuts power to the motor to prevent stripping.

Adaptable


Can be utilized across multiple anatomies with minimally invasive surgical options.

- 90-Degree collets (short & long) interface with the DePuy Synthes 90-Degree Screw driver for powered cordless insertion during a minimally invasive approach.




Designed to work with MatrixRIB[®] Fixation System, Titanium Sternal Fixation System, and MatrixMANDIBLE[™] Fixation System.

MatrixRIB[®] Fixation System



Titanium Sternal Fixation System



MatrixMANDIBLE[™] Plating System



MatrixPRO[®] II Driver

Adaptive Torque Limiting
Knows when to stop



Set 01.000.002

05.000.030	MatrixPRO II Driver Handpiece, qty of 2
05.000.031S	MatrixPRO II Driver Single Use Battery – Sterile, qty of 4
03.000.500	MatrixPRO II Driver Hex Collet, qty of 2
03.000.501	MatrixPRO II Driver J Latch Collet, qty of 1
03.000.502	MatrixPRO II Driver 90 Degree Collet, Short, qty of 1
03.000.503	MatrixPRO II Driver 90 Degree Collet, Long, qty of 1
60.000.015	MatrixPRO II Driver Graphic Case, qty of 1



For more information or to place an order, Please contact your DePuy Synthes sales consultant.

WARNING: In the USA, this product has labeling limitations. See package insert for complete information.

CAUTION: USA Law restricts these devices to sale by or on the order of a physician.

Not all products are currently available in all markets.



Distributed by:
(United States)
Synthes USA Products, LLC
1302 Wrights Lane East
West Chester, PA 19380
Telephone: (610) 719-5000
To order: (800) 523-0322
Fax: (610) 251-9056

Manufactured by:
(United States)
Pro-Dex, Inc.
2361 McGaw Avenue
Irvine, CA 92614
Telephone: (800) 562-6204
Fax: (866) 851-9339